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
I. TECHNICAL BASIS

The kit is based on a **Double antibody sandwich enzymatic immunoassay**. This technique shows big advantages from the rest of the techniques used till now: more sensitivity, more specificity, faster to do and easily automatizable.

Here we describe the technical basis of the assay: a rabbit polyclonal antibody against RHDV is coated on a solid support (polystyrene plate). When a sample containing

viral antigen is added, the viral particles will be caught by the antibody, After washing to remove all material not bound to the plate, we add a peroxidase conjugated antibody specific for RHDV. In the presence of viral antigen this conjugate will bind to it. We could demonstrate the presence of this labelled antibody by adding a chromogenic substrate that in presence of the enzyme will produce a colorimetric reaction.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive control and negative must be tested in a systematic way.
10. Substrate and stop solution must by handle whit care.

III. STORAGE OF COMPONENTS

Storage of the reagents: We recommend the storage of between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Macerate 1 gr of the sample (rabbit liver) in 2 ml of diluent. After a good homogenisation of the sample, centrifuge the homogenate for 10 min. at 1000 rpm. The supernatant obtained will be used as the assay sample.

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C.
- **Preparation of the conjugate to make immediately before use:**
Dilute the needed quantity of conjugate provided in the Kit 1/100 with diluent (e.g. 10µl of conjugated with 1 ml of diluent is more than enough quantity for a strip of 8 wells). Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.
- **Preparation of the control (+) and (-):**
Prepare a 1/100 dilution of the control samples provided with the kit in diluent (5 µl of control sample in 500 µl of diluent).

VII. TEST PROCEDURE

1. Take out the plate (or strips to be used) from the refrigerator and keep it 15 min at room temperature before starting the test.
2. Add 100 µl of positive control, diluted as it's specified in the previous instructions, to one well (e.i A1) and 100 µl of the negative control (prepared as is specified before) to another well (e.i. B1). Add 100 µl of each sample to test on each remainder wells. We recommend running the samples in duplicate wells.
Seal the plate and **incubate for 1 h at 37°C.**
3. Throw away the liquid contained in the plate on a recipient containing NaOH 0.1M. Wash 4 times following the procedure previously described.
4. Add 100 µl of conjugate (prepared as specified) to each well. Seal the plate and **incubate 1 hour at 25°C (room temperature).**
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate solution to each well. Keep the plate for **5 min at room temperature.**
7. Add 100 µl of stop solution to each well. It is recommended to add this reagent in the same order than was used for adding the substrate.
8. Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm.**

A. Validation of the Test:

The test is valid if:

$$OD_{450} \text{ Control (+)} > 1$$

$$OD_{450} \text{ Control (-)} < \text{cut off}$$

B. Results Interpretation:

The Cut Off value is the 15% of the OD₄₅₀ of the positive control sample:

$$\text{Cut of} = 0.15 \times OD_{450} \text{ Control (+)}$$

Samples with an OD higher than the Cut Off must be considered as positive to RHDV, and samples with OD value lower than the Cut Off must be considered as negative to RHDV. If samples in duplicate wells are being run, the OD value for the sample will be the average of OD obtained in both wells.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: