INGEZIM INFLUENZA A

Prod Ref: 1.0.FLU.K3

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I. TECHNICAL BASES OF THE KIT

The kit is based in the technique of blocking immunoenzymatic assay (Blocking ELISA) which is briefly described below:

Plates are coated with the antigen. After adding the sample to the well, if it contains specific antibodies against the virus, they will bind to the antigen adsorbed on plate while if the sample does not contain specific antibodies they will not bind to the antigen.

If we add a specific monoclonal antibody against the viral antigen coated to the plate (conjugated with

peroxidase), it will compete with the antibodies in the serum. If the sera samples contain specific antibodies, they will not allow the binding of the conjugate to the antigen; whereas if it does not contain specific antibodies, the Mab will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material, presence or absence of labelled Mab can be detected by adding the substrate (TMB) which, in presence of the peroxidase, will develop a colorimetric reaction.

II. PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. Avoid eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive and negative control must be tested in a systematic way.
10. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
11. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C. **Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES

Samples must be sera. To be used in the kit, they must be diluted ½ in diluent supplied in the kit. This dilution can be made directly in the well by adding 50µl of diluent and 50µl of sample.

VI. PREPARATION OF REAGENTS:

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionised water. When ready, this solution remains stable between +2°C and +8°C.
- **Positive and Negative controls:**
Controls must be tested as sample serum. Add 50µl of controls to the appropriate wells containing 50 µl of diluent.
¡ATTENTION! Once opened they are stable for one month between +2°C and +8°C. For long term storage, we recommend to make aliquots and store them at -20°C.

VII. TEST PROCEDURE

- All reagents except the conjugate must be brought to room temperature before use.
- If dilution of samples is going to be made directly in the well, add 50µl to the required wells and then, 50µl of controls and samples to the appropriated wells containing the diluent. Mix well contents to assure a perfect homogenization. Be careful to prevent cross-contamination. It is recommended to assay controls in duplicate. **Cover the plate with a lid and incubate at 37°C for 1 hour.**
- Shake out contents and wash 4 times as it is described previously.
- Add 100 µl of conjugate to each well. Cover the plate with a lid and incubate for **30 minutes at 37°C.**
- Shake out contents and wash 5 times as it is described previously.
- Add 100 µl of substrate to each well. The use of multichannel pipette is recommended to accurate time and uniformity. **Incubate at room temperature for 10 min.**
- Add 100 µl of stopping solution to each well in the same order and at the same speed as the substrate solution.
- Read the absorbance of each well using a ELISA reader with a filter of 450nm within 5 min. after the addition of stop solution.

VIII. READING AND RESULTS INTERPRETATION

A. Validation criteria:

The test is considered valid if:

$$\frac{\text{OD Negative Control}}{\text{OD Positive Control}} \geq 3$$

B. Results Interpretation:

There are two different ways for interpretation depending on the preferences:

- RESULTS EXPRESSED AS COMPETITION PERCENTAGE:**

$$\% \text{ Competition} = 100 - \frac{\text{OD Samples}}{\text{OD Negative Control}} \times 100$$

- Avian:**

- Samples will be considered positive if the percentage of competition is higher than 45 %.
- Samples will be considered negative if the percentage of competition is lower than 40%.
- In case of percentages between 40-45% a new evaluation is recommended with a new extraction after 2 weeks.

- Swine and Horse:**

- Samples will be considered positive if the percentage of competition is higher than 30 %.
- Samples will be considered negative if the percentage of competition is lower than 30%.
- In case of percentages between 25-35% a new evaluation is recommended with a new extraction after 2 weeks.

- RESULTS EXTRESED AS RATIO S/N**

$$S/N = \text{OD of the sample} / \text{OD of negative control}$$

- Avian:**

- Samples will be considered positive if the S/P < 0.55.
- Samples will be considered negative if the S/N > 0.60.
- In case of S/N between 0.55 y 0.60, a new evaluation is recommended with a new extraction after 2 weeks.

- Swine and Horse:**

- Samples will be considered positive if the S/N < 0.70.
- Samples will be considered negative if the S/N > 0.70.
- In case of S/P between 0.75 y 0.65, a new evaluation is recommended with a new extraction after 2 weeks.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: