INGEZIM ROTAVIRUS DAS

Prod Ref: 10.RT.K2

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I. THE TECHNICAL BASIS

he kit is based on a double antibody sandwich enzymatic immunoassay (**DAS-ELISA**). We make a brief description of the technique below:

We fix a monoclonal antibody (Mab) specific for Rotavirus on a solid support (Polystyrene plate). When a faecal sample contains the virus, the Mab will catch the viral antigen. After washing to eliminate all non-fixed material from the sample, we can detect

the presence of viral antigen using a rotavirus-specific Mab biotin-conjugated. After additional washing, we add a second conjugate (Streptavidin-peroxidase) which will allow the colour development in the presence of the substrate, during the last step of the assay.

In this way, the presence of colour means the presence of the virus in the sample and the absence of colour, the absence of viral antigen.

II. PRECAUTIONS AND WARNING FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. Distilled water must be used for reagents preparation.
9. For each utilisation of the kit, positive and negative control must be tested in a systematic way.
10. Substrate is very sensitive to light and contamination so never pipette directly from the substrate storage bottle. The necessary amounts for each assay must be always pour into a separate container.
11. Stop solution is a strong acid. Handle with care.

III. STORAGE OF THE KIT COMPONENTES:

Store all plates and reagents between +2°C and +8°C, except controls that must be stored at -20°C

IV. INFORMATION ABOUT THE WASHING STEPS:

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Homogenise 0.1 g of faecal sample in 1 ml of diluent (10% p/v).
Centrifuge at 1500 rpm for 10 minutes and collect the supernatant.

VI. PREPARATION OF REAGENTS:

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water. When ready this solution remains stable stored at 4°C.
- **Preparation of controls (+) and (-):**
Make 1/10 dilution in the supplied diluent before use.
- **Preparation of the conjugates 1 and 2: to make immediately before use.**
Dilute 1/100 with diluent:
 - The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.
 - The necessary and sufficient quantity of conjugate for an 8 wells strip is 10 µl of conjugate in 1 ml of diluent

VII. TEST PROCEDURE:

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of positive control, diluted as it's specified in the previous instructions, to two wells of the first row of the plate and 100 µl of the negative control to another two wells. Add 100µl of each sample to test, prepared following the previous instructions, on the remainder wells of the plate (is recommended to run samples in duplicate wells). Seal the plate and incubate for **60 min at 37°C**.
3. Wash 4 times following the procedure previously described.
4. Add 100 µl of conjugate 1, prepared following the previous instructions, to each well. Seal the plate and incubate **for 60 min at 37°C**
5. Wash 5 times following the described procedure.
6. Add 100 µl of conjugate 2, prepared following the previous instructions, to each well. Seal the plate and incubate **for 30 min at room temperature (20-25°C)**.
7. Wash 5 times following the described procedure
8. Add 100 µl of substrate. Keep the plate for **10 min at room temperature** in darkness.
9. Add 100 µl of stop solution to each well.
10. Read the absorbances of each well with an ELISA reader at 450 nm within 5 min after the addition of stop solution.

VIII. TEST PROCEDURE:

The reading must be done at 450 nm.

Validation of the plate:

For validation, OD values of the controls have

to fulfil these requirements:

OD of Positive control > 1.0
OD of Negative control < 0.2

Results Interpretation:

1º.- Calculate the corrected OD (cOD) for each sample:

$cOD = \text{sample OD} - \text{negative control OD}$

2º.- Interpretation :

Samples with a cOD higher than 0.3 are positives. If the cOD value is lower than 0.2, the sample is negative.

Samples with cOD ranged between these two values must be considered doubtful and should be retested.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: