INGEZIM GIARDIA CROM

Prod Ref: 10.GIA.K4

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

Esperar 10 minutos e interpretar:

I. TECHNICAL BASIS.

The technical basis of this test is a direct immunochromatography, where we are using a mixture of monoclonal antibodies specific to *Giardia lamblia*, which are able to detect all the shapes of the vital cycle of the parasites.

II. STORAGE OF THE REAGENTS

All the reagents must be kept between 4 and 25 °C.

III. PRECAUTIONS

1. Read carefully the instructions of use.
2. **Kept the reagents at room temperature before their use.**
3. Do not mix reagents or instruction from different kits.
4. Avoid any contamination of the reagents.
5. Do not use the kits after the expiry date.
6. There should be no eating, drinking or smoking where specimens or kit reagents are being handled.

IV. TEST PROCEDURE

1. Preparation of samples:

Unscrew the cap of the collection tube containing the diluent. Introduce the applicator stick into the faeces and collect a sample portion (30-50 mg). Reinsert the applicator stick into the tube and screw the cap. For liquid or semisolid stools, add 100 µl of stool using an appropriate pipette. Shake the collection tube to assure a right sample dilution.

2. Sample addition:

Open the aluminium bag just in the moment to realise the assay. Remove the test device from the bag and place it on a flat surface. Break the top part of the tube and dispense 3 drops into the window identified as "S" (sample). Wait 10 minutes and read.

3. Interpretation of results:


In order to be sure that the chromatographic assay has run properly, a line must appear in control area (window "C").

- POSITIVE RESULT: Lines in both areas "C" and "T".
- NEGATIVE RESULTS: Only line on the "C" area.
- NOT VALID RESULT: NO line on "C" area


Results interpreted after 10 minutes, may lead to the appearance of falsely positives.

The final diagnosis should not be based only on the result obtained with one test, it should be based on a correlation of tests results with typical clinical signs and symptoms.

Add 3 drops of diluted sample into the round window.

Wait 10 minutes and read the results.

Developed and Manufactured by:

Distributed by:

INMUNOLOGIA Y GENETICA APLICADA,
S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es