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I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the specific antigen on a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on plate. After washing to eliminate all non fixed material from the sera sample, we can detect

the presence of dog immunoglobulins using an specific peroxidase conjugate.

After addition of the substrate a colorimetric reaction will appear which could be measured by an spectrophotometer.

In this way the presence of colour means the presence of antibodies against the virus in the dog sera , and the absence of colour the absence of specific antibodies.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
9. Stop solution is a strong acid. Handle with care.
10. Avoid any contamination of the reagents of the Kit. Take special care of substrate because is very sensitive to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do

not maintain the plate on dry more time than strictly needed

- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

For the sera **TITRATION**, we recommend to make at least four two fold dilutions starting from the 1/50 dilution for each serum sample (1/50, 1/100, 1/200, 1/400).

VI. PREPARATION OF REAGENTS

WASHING SOLUTION:

Dilute one part of the concentrate washing solution provided in the kit into 9 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

PREPARATION OF THE CONJUGATE: TO MAKE IMMEDIATELY BEFORE USE.

Dilute 1/100 with diluent.

→ The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.

→ The necessary and sufficient quantity of conjugate for an 8 wells strip is 10µl of conjugate in 1 ml of diluent.

Prepare only the quantity needed for each time because the remainder volume has to be rejected.

Shake very well the solution prior to use.

CONTROL SERA.

Dilute as samples (1/50,1/100,1/200, 1/400).

VII. TEST PROCEDURE:

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of each positive control dilution on the first 4 wells of the first strip and 100 µl of the negative control dilutions on the remainder wells of the strip. Add 100 µl of each sera samples dilution to test, prepared following the previous instructions, on the remainder strips of the plate. Seal the plate and incubate for **10 minutes** at room temperature.
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and incubate for **10 minutes** at room temperature.
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate solution, to each well. Keep the plate for **5 min** at room temperature.
7. Add 100 µl of stop solution to each well.
8. Read the absorbances of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

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Results Interpretation:

The cut off value is equivalent to an O.D 450 nm = **0,4**

Positive samples: Samples with an OD higher than Cut off value.

Negative samples: Samples with an OD lower than the Cut off value.

For sera titration: The titre would be the maximum dilution with an O.D 450 nm higher than Cut off value.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: