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
I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on plate. After washing to eliminate all non-fixed material from the sera sample, we can detect the presence of cat immunoglobulins using a specific peroxidase conjugate. After addition of the substrate a colorimetric reaction will appear which could be measured by a spectrophotometer.

In this way the presence of colour means the presence of antibodies against the virus in the cat sera, and the absence of colour the absence of specific antibodies.

The aim of this kit is provide to the users with a reliable and automatizable diagnostic technique for this disease.

In our kit is important to remark the use of antigen obtained by the expression of the viral protein gp70 in baculovirus growth in insect cell cultures. This method warranties the total absence of infectivity in the kit.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
9. Stop solution is a strong acid. Handle with care.
10. Substrate is very sensitive to light and contamination. Do not introduce pipettes into the bottle. Pour off the volume exactly needed.

III. STORAGE:

Storage: All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

1. Throw out the content of the plate by a the possible mixture of the content from one well to another.
1. Dispense a volume of 300 µl of washing solution on each well.
1. Shake delicately the plate, avoiding the contamination between wells.
1. Turn over the plate brusquely to empty the wells.
1. Repeat the process as much times as is indicated on the instructions of the kit.
1. Prior to empty the content of the last brusque turn over of the plate to avoid
1. Do not maintain the plate on dry more time than strictly needed.
1. washing step, verify that the next.
1. After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

The sera samples need to be at 1/100 dilution in serum diluent to be tested. For titration we recommend to assay two fold dilutions from the first (1/200, 1/400, etc).

VI. PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the kit into 9 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

VII. TEST PROCEDURE:

1. Take out the plate (or strips) and the reagents and keep them at least 2 hours at room temperature before starting the test.
 2. Add 100 µl of positive and negative controls to duplicate wells. Add 100 µl of each dilution of sera samples to test, prepared following the Previous instructions, on the remainder wells of the plate (it is important to duplicate the controls and the samples tests)
Seal the plate and **incubate for 10 minutes at room temperature.**
 3. Wash 4 times following the procedure previously described.
 4. Add 100 µl of conjugate, to each well. Seal the plate and **incubate for 10 minutes at room temperature.**
 5. Wash 4 times following the described procedure.
 6. Add 100 µl of substrate solution, to each well. Keep the plate for **5 min at room temperature.**
 7. Add 100 µl of stop solution to each well.
 8. Read the absorbances of each well with an spectrophotometre at 450 nm
1. within 5 min after the addition of stop solution.

VIII. READING AND INTERPRETATION OF THE RESULTS:

The reading must be done with an ELISA reader at **450 nm**.

Kit Validation:

The test could be validated when:

-

O.D of positive control >1.

- O.D of negative control <0,25.

Results Interpretation:

Taking the negative control value, will be determined the following cut off:

-Negative cut off: O.D of negative control + 0,20

-Positive cut off: O.D of negative control + 0,25

Positive samples: Samples with an OD higher than the positive Cut Off value.

-Negative samples: Samples with a OD lower than the negative Cut Off value.

-Doubtful samples: Samples with an OD between both cut off. In these occasions, it is suggested to repeat the test after 3-4 weeks. If there is no seroconversion and its OD is not higher than the positive cut off in this second test, the sample will be considered negative.

If you run the samples in duplicate wells, the OD value for the sample will be the mean of both wells.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: