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## I. TECHNICAL BASIS

---

The kit is based on a **Double antibody sandwich enzymatic immunoassay**. This technique shows big advantages from the rest of the techniques used till now: more sensitivity, more specificity, faster to do and easily automatizable.

Here we describe the technical basis of the assay: a specific monoclonal antibody (Mab) against FeLV is coated on a solid support (polystyrene plate). When a sample containing viral antigen is added, the viral particles will be caught by the Mab. After we add a second specific Mab to FeLV.

This is a conjugate Mab (Mab-peroxidase). After incubation we could demonstrate the presence of this conjugate Mab by adding a chromogenic substrate that in presence of the enzyme will produce a colorimetric reaction.

The use of these very specific and stables Mabs, warrants the objectivity, specificity, safety and accuracy of the assay.

The coating Mab as well as the conjugate Mab, are very specific Mabs that allows the specific detection of FeLV antigen.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **iIMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. For each utilisation o f the kit, control positive and negative serum must be tested in a systematic way.
9. Stop solution is a strong acid. Handle with care.
10. Substrate is very sensitive to light an contamination. Do not introduce pipettes into the bottle.Pour off the volume extrictly needed.

## III. STORAGE

---

**Storage:** All reagents and plates must be stored between +2°C and +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

---

- ◆ The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well.
- ◆ After the incubation periods, the washing steps must be done following the next instructions:
- ◆ Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- ◆ Dispense a volume of 300 µl of washing solution on each well.
- ◆ Shake delicately the plate, avoiding the contamination between wells.
- ◆ Turn over the plate brusquely to empty the wells.
- ◆ Repeat the process as much times as is indicated on the instructions of the kit.
- ◆ Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use

- ◆ Do not maintain the plate on dry more time than strictly needed.
- ◆ After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

It is not necessary any previous dilution. Add directly 50 µl of the samples to each well.

## VI. PREPARATION OF REAGENTS

### Washing solution:

Dilute one part of the concentrate washing solution provided in the kit into 9 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

## VII. TEST PROCEDURE:

1. Take out the plate (or strips to be used) and keep it least 60 min. at room temperature before starting the test.
2. Add **50 µl** of positive control, and **50 µl** of the negative control in duplicate wells (prepared as specified before). Add **50 µl** of each sample to test on each remainder wells. We recommend running the samples in duplicate wells. Add **50 µl** of conjugate to each well (Controls and samples).
3. Wash 4 times following the described procedure.
4. Add 100 µl of substrate solution, to each well. Keep the plate for **5 min at room temperature**.
5. Add 100 µl of stop solution to each well.
6. Read the absorbances of each well with an ELISA reader at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND INTERPRETATION OF THE RESULTS:

The reading must be done with an ELISA reader at **450 nm**.

### Kit Validation:

The test could be validated when:

- O.D of positive control >1.**
- O.D of negative control <0,20.**

### Results Interpretation:

Taking the negative control value, will be determined the following cut off:

- Negative cut off: O.D of negative control + 0,20**
- Positive cut off: O.D of negative control + 0,25**

### Positive samples:

Samples with an OD higher than the positive Cut Off value.

### Negative samples:

Samples with a OD lower than the negative Cut Off value.

### Doubtful samples:

Samples with an OD between both cut off. In these occasions, it is suggested to repeat the test after 3-4 weeks. If OD is not higher than the positive cut off in this second test, the sample will be considered negative.

If you run the samples in duplicate wells, the OD value for the sample will be the mean of both wells.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in \_\_\_\_\_ by: