INGEZIM FeLV VET

Prod Ref: 16.FLV.K9-32

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I. TECHNICAL BASIS

The kit is based on a **Double antibody sandwich enzymatic immunoassay**. This technique shows big advantages from the rest of the technique used till now: more sensitivity, more specificity, faster to do and easily automatizable.

Here we describe the technical basis of the assay: and specific monoclonal antibody (Mab) against FeLV is coated on a solid support (polystyrene plate). When a sample containing viral antigen is added, the Mab will catch the viral particles. After, we add a second specific Mab to FeLV that is

conjugated to HRPO (Mab-peroxidase). After incubation we could demonstrate the presence of this conjugate Mab by adding a chromogenic substrate that in presence of the enzyme will produce a colorimetric reaction.

The use of these very specific and stable Mabs, warrants the objectivity, specificity, safety and accuracy of the assay.

The coating Mab as well as the conjugate Mab, are very specific Mabs that allows the specific detection of FeLV.

II. PRECAUTIONS AND WARNING FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **(IMPORTANT)**.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
8. Use a new pipette for each serum sample.

III. STORAGE CONDITIONS:

All reagents and components, must be stored at 4 °C.

IV. INFORMATION ABOUT THE WASHING STEPS:

The washing steps could be done using a washing solution bottle or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instruction:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense the washing solution (prepared following instructions), proportionally to the wells used.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the kit.

- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
- Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

It is not necessary any previous dilution. Add directly 1 drop of the samples to each well.

VI. PREPARATION OF REAGENTS:

All the reagents are ready to use, without any previous manipulation. Only the washing solution is 10 x concentrated.

WASHING SOLUTION PREPARATION: ®
Dilute one part of the concentrate washing solution provided in the kit into 9 parts of

distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

Example: 60 ml. of concentrated washing solution into 540 ml. of distilled or deionized water.

VII. TEST PROCEDURE:

1. Take out all the reagents and the strips or wells to be used and keep it at room temperature before starting the test.
2. Add **1 drop** of positive and negative control on two different wells. Add **1 drop** of each sera sample to test on each remainder wells using one different pipette for each one.
Add **1 drop** of conjugate to each well (Controls and samples). shake the frame softly for some seconds. Seal the plate and incubate for **5 minutes at room temperature.**
3. **Wash 5 or more times** following the described procedure.
4. Add **1 drop** of chromogen solution, to each well and **1 drop** of substrate solution in the same way. Keep the plate for **5 min at room temperature.**
Read the results.

VIII. READING AND INTERPRETAION OF THE RESULTS:

The reading must be done by visual colour turn detection.

Kit validation:

We take the controls like reference.

- **Negative control:** There is no colour turn detection (transparent).
- **Positive control:** There is colour turn detection (blue).

Results Interpretation:

- **Negative samples:** There is no colour turn detection.
- **Positive samples:** There is colour turn detection.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: