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INGEZIM RABBIT 17.RHD.K1

I. TECHNICAL BASIS

The kit is based on the indirect immunoenzymatic assay technique described below:

The antigen is fixed in a solid support (polystyrene plate). The antigen is recombinant protein from RHDV obtained using the baculovirus system of expression.

The serum sample is added to each well. After incubation, a peroxidase conjugate (Protein A) is added. If serum sample contains antibodies specific of RHDV, the conjugate will bind to them, whereas if it does not contain specific antibodies no binding will be appreciated.

After washing the plate to eliminate all non fixed material, we can detect presence or absence of

conjugate by adding the substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

The assay has been compared with hemagglutination inhibition technique with the highest ratios of sensitivity and specificity.

The kit could be used both for screening and for titration purposes. We have adapted the conditions for verifying in a very ease way the level of protection of a vaccinated herd. The process is based on experimental infection trials. Using this procedure the verification of the protection level could be done using the screening scheme at a single dilution.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
10. Substrate and substrate buffer must by handle whit care.

III. STORAGE OF COMPONENTS

Store plates and reagents between +2°C and +8°C except the conjugate that must be placed at -20°C as soon as the kit is received. The conjugate is supplied in a glycerol solution that warrants the conservation of the reagent in liquid phase even under -20°C.

ATTENTION! Stop solution may precipitate at this temperature, if it happens, warm at 37°C till total dissolution

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

FOR TITRATION PURPOSES:

We recommend to make at least four 2 fold dilutions from 1/200 (1/200, 1/400, 1/800 and 1/1600), but depending on the goal, higher dilutions could be done. Dilution must be done in the supplied Diluent.

FOR SCREENING PRUPOSES:

Before assaying, sera samples must be diluted 1/200 in the supplied Diluent (5 µl of serum sample to 1 ml with diluent).

This is also the procedure recommended for the verification of the protection of animals..

VI. PREPARATION OF REAGENTS

- Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water(40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C.
- Diluent:**
Dilute one part of the concentrate diluent provided in the Kit with 4 parts of distilled or deionized water.
- Control Sera:**
Controls are ready to use.
- Preparation of the conjugate (to make immediately before use):**
Dilute the needed quantity of conjugate provided in the Kit 1/100 with the provided diluent:
 - The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of Conjugate with 11 ml of diluent.
 - The necessary and sufficient quantity of Conjugate for an eight wells strip is 10 µl of conjugate with 1 ml of conjugate diluent. Shake very well the solution before use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.**Substrate:** Substrate is ready to use.

VII. TEST PROCEDURE

- Take out the plate (or strips to be used) from the refrigerator and keep it 15 min at room temperature before starting the test.
- Add 100 µl of positive control serum and 100 µl of negative control serum to A1 and B1 wells respectively.
Add 100 µl of each sera sample (prepared as specified) to each remainder wells of the assay plate. For high accuracy it is recommended to run samples and controls in duplicate. Seal the plate and **incubate for 1 hour at 37°C**.
- Wash 3 times following previous instructions
- Add 100 µl of conjugate (prepared as specified) to each well. Seal the plate and **incubate 1 hour at 25°C (room temperature)**.
- Wash 5 times following the described procedure.
- Add 100 µl of substrate, to each well. Keep the plate for **10 min at room temperature**.
- Add 100 µl of stop solution to each well. It is recommended to add this reagent in the same order than was used for adding the substrate.
- Read the OD of each well at 405 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **405 nm**.

A. Validation of the Test:

To validate the assay, the ratio OD (Positive Control) / OD (Negative control) must be higher than 10.

B. Interpretation of the Results:

A. Cut off calculation:

Cut off Positive/Negative is 0,300.

Hence, samples with an OD higher than 0,300 must be considered as positives and samples with OD lower than 0,300 must be considered as negatives.

- SCREENING PURPOSES:**

All samples with OD higher than 0,300 must be considered as positives, and samples with OD lower than this value, must be considered as negatives.

- TITRATION PURPOSES:**

The sera titer will be the highest dilution with an OD higher than 0.3.

- PROTECTION VERIFICATION:**

After running correlation trials for titer / OD value, we can estimate than samples with an OD value (at a 1/200 dilution) higher than 0.9, are enough protected to experimental infection with RHDV. Anyhow this is a general orientation that could vary depending on the vaccination plans, etc. In order to establish your particular level of protection, we recommend running your own seroprofiles using the titration procedure

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by:

IT-73840
IT-73780

ISO 14001:2015
9191.INGE

ISO 9001:2015
9175.ING2