INGEZIM ASFV CROM Ag

Prod Ref: 11.ASF.K42

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

NEGATIVO

I. TECHNICAL BASIS

Ingezim® ASFVCROM Ag is an immunochromatographic assay for the detection of ASFV in blood samples. The test is based on the use of two different coloured latex microspheres: black microspheres that are coated with a specific MAb to ASFV, and blue microspheres, which are used as test control.

On the membrane, two lines are printed: the test line (T) has a specific MAb to ASFV and the control line (C) has a specific MAb for the control protein.

If the case of a positive sample, the virus binds to the black beads conjugated to an anti-ASFV MAb. The immune complex then migrates through the membrane by capillarity and is captured again by the anti-AFSV MAb absorbed on the test line, resulting in the appearance of a black line. The presence of the control line serves as validity of the test, indicating that the immunochromatography has been performed correctly.

I. STORAGE

All the reagents must be kept between +4°C and +25 °C.

1. **PRECAUTIONS** Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration date and do not mix components from different batches.
6. Do not eat, drink, or smoke where specimens or kit reagents are being handle.

II. SAMPLES

The test has been designed for use with blood samples. Blood samples must be fresh or refrigerated up to 4 days at 2-8 °C, and collected with any anticoagulant (EDTA, heparin, etc.). Blood samples without any anticoagulant, could contain micro-clots that can block the device and result in nonspecific reactions.

III. TEST PROCEDURE

- **Allow all kit components and samples to reach room temperature, before use.**
- Remove the device from the aluminium pouch and discard the desiccant.
- Add 20 µl of the whole blood (with anticoagulant) to the round window, and wait for one minute.
- Add 3-4 drops, **drop by drop and waiting for the total absorption between drops**, of the running buffer on the same round window. Wait 10 minutes to readout the results.

IT IS VERY IMPORTANT TO FOLLOW THESE INSTRUCTIONS AND RESPECT THE INDICATED VOLUMES TO OBTAIN VALID RESULTS

IV. INTERPRETATION OF RESULTS


Interpretation of results should be done as follows, 10 minutes after the addition of the sample:

NEGATIVE: A single blue line appears in the result area.


POSITIVE: A blue line and a second black line appear in the result area. Colour intensity may vary depending on the antigen concentration.

INVALID: If no blue line appears, the test will be invalid. Repeat the test using a new device.


- **No results must be considered after 10 minute.**


1. Add 20 μL of sample to the round window

2. Add 3-4 drops of running buffer

POSITIVE

NEGATIVE

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENÉTICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91 368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: