INGEZIM AHSV CROM ANTICUERPO

Prod Ref: 14.AHS.K41

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

**POSITIVO****NEGATIVO**

I. TECHNICAL BASIS

The test is based on the technique of immunochromatography, a migration technique that uses purified VP7 protein of AHSV, and is capable of detecting specific antibodies against AHSV in porcine serum samples.

The diagnostic device consists of two windows:

- Sample window: it contains VP7 protein of AHSV and a protein control, coated to coloured latex particles.
- Result window: it contains a test line (T) formed by the VP7 protein of AHSV and a control line (C) formed by a MAb specific for the protein control.

If the sample contains antibodies against the VP7 protein of AHSV, they will react with the red particles conjugated to VP7 protein. The latex-protein-antibody complex migrates through the membrane and reacts again with the VP7 protein adsorbed in the test line resulting in the appearance of a red/pink. The appearance of a blue line in the control area (C) indicates that the chromatography has been correctly performed.

II. STORAGE:

All the reagents must be kept between 4 and 25 °C.

III. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle

IV. TEST PROCEDURE


Reagents and devices must be at room temperature when used for testing. Open the aluminium bag just in the moment to realise the assay. Remove the test device from the bag and place it on a flat surface.

Dispense 15 µl of the serum sample into the sample well (round window).


Shake vigorously the dropper bottle with running buffer and add 6-7 drops into the sample window (**wait total liquid absorption into the membrane between drop and drop**).

Read results at 10 minutes.

INGEZIM AHSV CROM *antibody*
14.AHS.K41

1. Add 15 μ l of serum


2. Add 6-7 drops of diluent buffer

V. INTERPRETATION OF RESULTS


In order to be sure that the chromatographic assay has run properly, a blue line must appear in control area (C).

- If appears a red/pink line in the test area (T) the sample assayed must be reported as POSITIVE to AHSV.
- If there is no line in the test area (T) and there is a blue line in the control area (C) the sample must be reported as negative to AHSV.
- NO VALID TEST: If there is no blue line into the control area (C) the assay can't be considered valid.

No result must be considered after 10 minutes

POSITIVE

NEGATIVE

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by:

IT-73840
IT-73780

ISO 14001:2015
9191.INGE

ISO 9001:2015
9175.ING2