INGEZIM FVR IgM

Prod Ref: 13.FVR.K2

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I. TECHNICAL BASIS

This kit is based on a capture enzymatic immunoassay. The technique is briefly described below:

A Monoclonal antibody (MAb) specific to ruminant IgMs is fixed on a solid support (polystyrene plate). The IgMs present in the samples will be captured by the MAb adsorbed on the plate. After a washing step to eliminate all non-fixed material from the serum sample, a viral antigen of the Rift Valley Virus (recombinant N protein), peroxidase conjugated, is added. If the serum sample contains specific

immunoglobulins M to FVR this viral conjugate will be captured.

After a washing step and adding the substrate a colorimetric reaction will develop, which can be measured by an ELISA reader.

The development of colour means that there are IgM antibodies against the virus present in the ruminant sera, whilst the absence of colour indicates the absence of these specific antibodies.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Do not use components after expiry dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking while handling reagents or samples.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. Systematically include a positive control and a negative serum each time the assay is run.
9. The stop solution is a strong acid and therefore must also be handled with care.
10. The substrate must be handled with care, as it is very sensitive to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2 °C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µl of washing solution in each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

V. PREPARATION OF SAMPLES:

Sera samples must be tested at 1/5 dilution in the serum diluent. This dilution can be made directly on the plate by adding 80µl of the diluent and 20µl of the sample to each well. Mix gently to homogenize the solution.

VI. PREPARATION OF REAGENTS

Washing solution 25x:

Dilute one part of the concentrated washing solution provided in the kit into 24 parts of distilled or deionized water. (40ml of the concentrated solution in 960ml of distilled or deionised water). Once prepared, this solution remains stable when stored between +2°C and +8°C.

Preparation of the sera diluent

The diluent is ready to use. **Do not dilute**

Preparation of the control sera:

Controls are ready to use. Do not dilute.

Preparation of the viral antigen peroxidase conjugated:

The antigen is ready to use. Do not dilute.

VII. TEST PROCEDURE:

- Prior to starting the test, bring all reagents to room temperature (20-25°C) (except for the conjugate, which must be brought to RT right before its use).
- Incubation of serum samples:
 - Add 80 µl of the diluent to each well and then add 20µl of the serum to each well.
 - Homogenize the content of the wells by gently shaking the plate.
 - Add 100µl of the ready to use controls.
 - Seal the plate and **incubate for 1h at room temperature (20-25°C)**.
- Wash 4 following the procedure previously described.
- Add 100 µl of the conjugate to each well. Seal the plate and incubate for **30 min at room temperature (20-25°C)**.
- Wash 5 times following the procedure previously described.
- Add 100 µl of the substrate solution, to each well. Keep the plate for **15 min at room temperature**.
- Add 100 µl of the stop solution to each well.
- Read the OD of each well at 450 nm within the following 5 min after the addition of the stop solution.

VIII. READING AND RESULT INTERPRETATION

Validation of the test:

The test is considered valid when the OD of the positive control at 450 nm minus the OD of the negative control is higher than 0.25 and the OD of the negative control at 450 nm is lower than 0.25.

Interpretation of the results:

The S/P must be calculated as follows:

Samples with an S/P higher than or equal to 0.4 should be considered as **positive** to FVR IgM antibodies.

Samples with an S/P lower than 0.4 should be considered as **negative** to FVR IgM antibodies.

(Sample OD)

S/P = _____

(Positive Control OD)

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by:

IT-73840
IT-73780
ISO 14001:2015
9191.INGE
ISO 9001:2015
9175.ING2