## INGEZIM Schmallenberg COMPAC 2.0

Prod Ref: 13.SBV.K3

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## INGEZIM SCHMALLENBERG Compac 2.0 13.SBV.K3

### I. TECHNICAL BASIS

The Kit has been designed to detect, in a very easy way, specific antibodies to SBV in sera samples. The assay has been developed to use mainly in sheep but it can be used for other ruminants. This Kit is based on a blocking enzymatic immunoassay. We make a brief description of the technique bellow:

The SBV antigen is fixed in a solid support (polystyrene plate). The serum sample is added to each well. After incubation, we add a specific monoclonal antibody (peroxidase conjugated) specific to SBV. If the serum sample do not contains antibodies against the virus, the conjugate will bind to

the plate, whereas if it contains specific antibodies to SBV, these antibodies will bind the antigen on the plate blocking the binding of the conjugate, that's why no binding of the conjugate will be appreciated.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

Due to the high specificity of the monoclonal used as conjugate, the assay has a high specificity and sensitivity index.

### II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handle whit care, it is very sensible to light and contamination.

### III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

### IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

### V. PREPARATION OF SAMPLES:

A 1/10 dilution must be made before assaying the sample. This dilution could be made directly on each

well by adding 90µl of serum diluent and 10µl of serum sample and mixing properly.

## VI. PREPARATION OF REAGENTS

### • *Washing solution:*

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml de concentrated solution and 960 ml of water). Once ready, this solution remains stable between +2°C and +8°C.

### • *Preparation of controls:*

Controls are ready to use (**DO NOT DILUTE**) and must be handled adding 100 µl/well.

### • *Preparation of the conjugate:*

The conjugate is ready to use (**DO NOT DILUTE**) and must be handled adding 100 µl/well.

## VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Add 90 µl of diluent to all the wells to be used with samples.
  - Add 10 µl of serum samples to each well plate. Shake carefully the plate to help the correct homogenisation.
  - Add 100 µl of Positive and Negative controls in duplicated wells.

Seal the plate and **incubate for 1 hour at room temperature (22-25°C)**.
- Wash 5 times following the described procedure
- Add 100 µl of conjugate to each well. Shake carefully the plate to help the correct homogenisation of the components mixed on each well. Seal the plate and **incubate for 30 min at room temperature (22-25°C)**
- Wash 5 times following the described procedure.
- Add 100 µl of substrate solution, to each well. Keep the plate for **15 min at room temperature, in a dark place.**
- Add 100 µl/well of stop solution. NOTE: Stop solution must be added in the order in which the substrate was added.
- Read the OD of each well at 450 nm with a spectrophotometer within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading of results must be done at **450 nm** with a spectrophotometer.

### A. Validation criteria:

- OD of Negative control >0.750
- To validate the assay, the ratio OD (positive Control) / OD (negative control) must be lower than 0.3

### B. Results Interpretation:

When samples were run in duplicate, arithmetic mean of both values must be calculated for each sample.

Calculate the blocking percentage for each sample as follows:

$$\%B = (1 - (DO_{\text{Muestra}} / DO_{\text{Control Neg}})) * 100$$

- All samples with blocking percentage lower than or equal to 50% must be considered Negative.
- All samples with blocking percentage higher than or equal to 55% will be considered Positive.
- Samples with blocking percentage between both should be considered doubtful. For these samples a new assay is recommended to be run after 3 weeks.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in

by:

A large empty rectangular box with a blue border, intended for the distributor's name.