

## INgezim Schmallenberg Compac

R.13.SBV.K3

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

## INgezim Schmallenberg Compac

R.13.SBV.K3

**INgezim SBV COMPAC** is a multispecies enzymatic immunoassay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to the N protein of Schmallenberg virus (SBV).

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with SBV antigen. Serum samples are added and incubated.
2. If the samples contain specific antibodies to SBV N protein they will bind to the antigen.
3. When a MAb-PO specific of N protein of SBV is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection of specific antibodies to SBV N protein in domestic and wild ruminants serum samples.

### INTERPRETATION OF THE RESULTS

A Cut off regarding an **inhibition percentage (IP)** is established. Samples will be considered **Positive** if their IP value is equal to or higher than the positive cut off. Samples showing an IP within the Cut offs will be considered doubtful.

## VALIDACIÓN

### 1. Experimental animals:

597 ovine sera (48 vaccinated positive to seroneutralization (SN), 3 infected positive to SN, 546 free of SBV), 92 goat sera free of SBV, 132 bovine sera (40 vaccinated positive to SN, 92 free of SBV).

**Sensitivity respect to SN was 96% and specificity higher than 99%**

### 2. Field animals:

186 mouflons, 109 wild boars, 299 fallows, 49 roe deer, 1053 deer, 31 mountain goats, 213 cattle, 89 goats, 203 sheep.

**Sensitivity of the assay respect to SN was higher than 99% and specificity within the range of 95% and 100% depending on the specie studied.**

**When general calculations are made (experimental and field animals), results obtained depending on the source of the animals (domestic or wild) are shown in the table below:**

distribution	sensitivity (%)	specificity (%)
Domestic	99	99
Wild	>99,5	99

### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

PRODUCT MANUFACTURED BY INGENASA  
REGISTRATION NUMBER: 3286 RD

SHELF LIFE: **12 months**  
Stored at 2°C-8°C

Ed.050619