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
I. TECHNICAL BASIS

This kit is based on a capture enzymatic immunoassay. We make a brief description of the technique below:

A Monoclonal antibody (MAb) specific of equine IgM is fixed on a solid support (polystyrene plate). Samples are placed in two wells of microplates. If sample contains IgM antibodies they will be captured by the MAb adsorbed on the plate. After washing to eliminate all non fixed material from serum sample, viral antigen of West Nile Virus (WNVAP) should be added to 1 well and Control Negative Antigen (CNA)

to the second well .If serum sample contains specific immunoglobulins M to WNV they will be captured.

After additional washing peroxidase-conjugated Mab specific to WNV E protein (domain III) should be added. After washing and addition of the substrate a colorimetric reaction will be developed which can be measured by an ELISA reader.

In this way presence of colour means presence of antibodies IgM against the virus in equine sera, and absence of colour, absence of specific IgM antibodies.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. Do not eat, drink or smoke where specimens or kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. The positive and negative control supplied with the kit must be added to each series of analyses.
9. Stop solution is a strong acid. Handle with care.
10. Substrate is very sensitive to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2 °C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do

not maintain the plate on dry more time than strictly needed.

- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Samples must be tested in two wells. Enough quantity for two wells should be prepared. Sera samples must be tested at dilution 1/5 in serum diluent. This dilution can be made on the assay plate directly by adding 40µl of diluent and 10µl of sample to each well. Mix gently to homogenate the solution.

VI. PREPARATION OF REAGENTS

Washing solution 25x:

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water. (40ml of the concentrated solution in 960ml of distilled or deionised water) When ready this solution remains stable when is stored between +2°C and +8°C.

Preparation of the sera diluent

Diluent is ready to use. **Do not dilute**

Preparation of the control sera:

Controls are ready to use. **Do not dilute.**

Preparation of the viral antigen and the control antigen:

Antigens are ready to use. **Do not dilute.**

Preparation of the conjugate

Conjugate is ready to use. **Do not dilute.**

VII. TEST PROCEDURE:

1. It is recommended to bring the reagents to be used to RT (20-25°C). Let the conjugate reach the RT just before be used.

2. Incubation of serum samples:

ALL THE SAMPLES MUST BE RUN BY DUPLICATE. (see Picture 1)

- Add 40 µl of diluent to each well and then add 10µl of the serum to each well (two wells per sample).
 - Homogenise content of wells by gentle shaking of the plate.
 - Add 50µl of controls ready to use (two wells per control).
 - Seal the plate and **incubate for 1h at room temperature (20-25°C) in a moist chamber**
3. Wash 3 times following the previously described procedure.

4. Dispense 50µl of WNV positive antigen (APWNV red colour) in the well containing the first duplicate, and 50µl of negative antigen control (ANC green colour) in the well containing the second duplicate. (see figure 1)

5. Seal the plate and **incubate for 16-24h at room temperature (20-25°C) in a moist chamber**

6. Wash 5 times following the described procedure

7. Add 50 µl of conjugate to each well. Seal the plate and incubate for **1h at room temperature (20-25°C) in a moist chamber.**

8. Wash 5 times following the described procedure.

9. Add 50 µl of substrate solution, to each well. Keep the plate for 15 min at room temperature.

10. Add 50 µl of stop solution to each well.

11. Read the OD with an ELISA reader at 450 nm (within 5 min after stopped)

1.STEP Run the samples by duplicate
2*STEP Add positive antigen and negative antigen to each sample

APWN: Positive Antigen to WNV
ANC: Negative Antigen control

CN: Negative Control
CP: Positive Control

M1: Sample 1 by duplicate in wells C1 & C2
M2: Sample 2 by duplicate in wells D1 & D2

Picture 1: ELISA Scheme

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

The OD obtained with negative antigen must be subtracted from the OD obtained with positive antigen (sample and controls)

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B. Results Interpretation:

- *Positive and negative samples*

Sample should be considered positive for IgM antibodies to WNV when

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Sample should be considered negative for IgM antibodies to WNV when

ODcm: (OD POS Ag – OD NEG Ag) ≤ 0.30

- *Doubtful samples*

Samples showing values between 0.30 and 0.35 ("Grey zone") should be considered doubtful and retested. If the result is the same we recommend to assay a new extraction taken off 15 days after the first one.

All positive and doubtful (grey zone) samples must be confirmed by seroneutralization.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: