ASFV

Anticuerpo Monoclonal específico de vp72 de VPPA VP72 ASFV specific Monoclonal Antibody

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The ASF is one of the most complex viral diseases that can affect swine. The porcine species are the only species found to be naturally susceptible to the disease. It produces clinical signs ranging from an acute form to a subacute and/or chronic form characterized by high fever, high mortality rates (in acute form) extensive hemorrhages in internal organs, and functional alterations of the respiratory and nervous systems.

Clinically, ASF may resemble a variety of other swine hemorrhagic diseases, especially Classical Fever. A laboratory test is required to establish a correct diagnosis.

ASF is endemic in many countries of sub-Saharan Africa. In European countries, where ASF is still present, wild pigs (boars, etc.) play a major role in the epidemiology of the disease.

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
The hybridoma which produces the monoclonal antibody has been obtained by the fusion of lymphocytes from Balb/c mice's spleen with myeloma X63/Ag8653 cells. The IgG has been purified by HPLC, showing a purity of 99%. Its isotype is IgG₁ and it is specific to ASFV VP72.

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- Seed 40.000 Vero cell / well in slides.
- The following day, remove the culture medium and infect with Lisbon 60 strain of ASFV.
- Incubate 2 h. at 37°C and remove the inoculum.
- Add culture medium
- Incubate 24 h. and remove the culture medium.
- Wash once with PBS.
- Dry completely at R.T.
- Saturate with PBS: BSA 1%
- Add first antibody and incubate in wet chamber at R.T. (Make dilutions in PBS: BSA 0.25%) (See recommended use conditions for dilution. Note that the supernatant is supplied ready to use)-
- Wash by immersion.
- Add the second antibody labelled with fluorescein diluted in PBS: BSA 0.25%.
- Incubate in wet chamber at R.T.
- Wash by immersion.
- Prepare the sample with glycerol:PBS (9:1)
- Fix the slide.

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